



Zfx™ GenTek™ Abutment Blanks

Kompatible Abutment Halter

Zfx™ Abutment Halter	Artikel-Nr.
Zfx™ Pre-Abutment Blank Halter, Standard ¹	ZFX02PMSH
Zfx™ Pre-Abutment Blank Halter, Multi-Blank für Zfx™ Inhouse5x	ZFX02002069
Zfx™ Pre-Abutment Blank Adapter für Glaskeramik Halter ²	ZFX02PMA

Medentika® PreFace® Abutment Halter ³	Artikel-Nr.
PreFace® Abutment Halter Datron D5	
PreFace® Abutment Halter Dental Concept Systems DC1, DC5	
PreFace® Abutment Halter Dental Concept Systems DC5	
PreFace® Abutment Halter imes-icore 140i	
PreFace® Abutment Halter imes-icore 350i, 450i, 550i, 650i, 850i	
PreFace® Abutment Halter imes-icore 750i	
PreFace® Abutment Halter MB Maschinen Cobra Mill 5A1, 5M	
PreFace® Abutment Halter Primacon PFM 24 mediMill	
PreFace® Abutment Halter R + K Organical Multi5X, 5XT-M, Desktop 8	
PreFace® Abutment Halter Rödgers RXD	
PreFace® Abutment Halter Sirona inLab MC X5	
PreFace® Abutment Halter VHF N4	
PreFace® Abutment Halter VHF S2	
PreFace® Abutment Halter Wissner Gamma 202	

PreFace® Abutment Halter müssen direkt beim Maschinenhersteller bestellt werden.

¹ Der Pre-Abutment Halter, Standard kann auf allen Maschinen verwendet werden, für die ein Standard Blank-Halter mit einem Durchmesser von 98 mm und einer Höhe von 10 mm verwendet werden kann.

² Der Zfx™ Pre-Abutment Adapter kann auf allen Spannsystemen verwendet werden, die eine 10 x 10 mm Verbindung einspannen können.

³ Zfx™ GenTek™ Pre Abutment Blanks können mit jeder Fräsmaschine gefräst werden, die mit einem Medentika® Halter ausgestattet ist. Diese Liste ist keine vollständige Liste kompatibler Fräsmaschinen.