




Product Catalog  
2023



# GenTek™

## Genuine Restorative Components





# GenTek Open Connectivity

For even more flexibility.

## Genuine

GenTek Pre-Abutment Blanks are compatible with Medentika® holders (CAM update required).



## NEW

### GenTek TiBases


Achieve efficiency to your CAD/CAM restoration workflow with a singular TiBase that provides flexibility of channel angulation and adjustable post height options.

exocad

RealGUIDE™  
UNIVERSAL OPEN SYSTEM

dental wings

3shape




## Free

Libraries for the GenTek restorative components are available for Zfx CAD, RealGUIDE™, 3Shape, exocad, and Dental Wings CAD software.

Free download at [zfx-dental.com](http://zfx-dental.com)



SCAN & INFORM  
More at  
[www.zfx-dental.com](http://www.zfx-dental.com)







## Table of Contents

<b>NEW: GenTek TiBases</b> .....	<b>4</b>
<b>GenTek Pre-Abutment Blanks</b> .....	<b>8</b>
<b>GenTek TiBases – Milling Block Compatible</b> .....	<b>12</b>
<b>GenTek TiBases – ZimVie Multi-Unit Abutment Compatible</b> .....	<b>20</b>
<b>GenTek Scanbodies</b> .....	<b>22</b>
<b>GenTek Digital Scan Analogs</b> .....	<b>28</b>

## Legend

**∅ Restorative Platform (mm)**  
Color code for Eztec, TSX and Tapered Screw-Vent





	2.9 mmD
	3.5 mmD
	4.5 mmD
	5.7 mmD

## Placement Tool for GenTek Digital Scan Analogs Compatibility Guide

Placement Tool	Laboratory Analog
M 14	Low Profile Abutment Lab Analogs
M 16	Certain and Eztec Lab Analogs
M 18	TSV and Tapered Abutment Lab Analogs
M 20	External Hex Lab Analogs



**∅ Restorative Platform (mm)**  
Color code for Certain and External Hex

	3.4 mm
	4.0–6.0 mm
	5.0 mm
	6.0 mm

## GenTek Calibration Tool for Milling Machines

Description	Part No.
GenTek Calibration Abutment Blank	ZFX18-CALIB



## GenTek Pre-Abutment Blank Holder

Description	Part No.
Zfx Pre-Abutment Blank Holder, Standard	ZFX02PMSH
Zfx Pre-Abutment Blank Holder, Multi-Blank for Zfx Inhouse5x	ZFX02002069
Zfx Pre-Abutment Blank Adaptor for Glass Ceramic Holder	ZFX02PMA



## Milling Block Compatibility Guide<sup>1</sup>

ZimVie Restorative Platform Connection	Block Interface Size
Tapered Screw-Vent and Trabecular Metal	L
Eztec	S
Certain	S / L <sup>2</sup>
External Hex	L



<sup>1</sup> Not applicable for the GenTek Angulated Screw Channel TiBase

<sup>2</sup> Restorative Platform 3.4 mm is compatible with interface size S  
Restorative Platform 4.1 and 5.0 mm are compatible with interface size L

# NEW: GenTek TiBases

## ZimVie Dental Certain Restorative Platform Connection

### TiBase Engaging

*incl. Gold-Tite Hexalobular Screw and Gluing Screw*



Post Height <sup>1</sup> (mm)	Restorative Platform (mm)	Collar Height (mm)	Part No.	
7	3.4	0.3	ZFX11ZB-CE34ASZ03E	
		1.3	ZFX11ZB-CE34AS13E	
		2.6	ZFX11ZB-CE34AS26E	
	4.1	4.1	0.3	ZFX11ZB-CE41ASZ03E
			1.3	ZFX11ZB-CE41AS13E
			2.6	ZFX11ZB-CE41AS26E
	5.0	5.0	0.3	ZFX11ZB-CE50ASZ03E
			1.3	ZFX11ZB-CE50AS13E
			2.6	ZFX11ZB-CE50AS26E

### Gold-Tite® Hexalobular Screw

### Hexalobular Try-In Screw



Part No.		Part No.	
ZFX09-ZB-CE-HLGTE	20 Ncm	ZFX09-ZB-CE-HLTISE	Hand-tighten

### TiBase Non-Engaging

*incl. Gold-Tite Hexalobular Screw and Gluing Screw*



Post Height <sup>1</sup> (mm)	Restorative Platform (mm)	Collar Height (mm)	Part No.
7	3.4	0.3	ZFX11ZB-CE34ASZ03N
		1.3	ZFX11ZB-CE34AS13N
	4.1	0.3	ZFX11ZB-CE41ASZ03N
		1.3	ZFX11ZB-CE41AS13N
	5.0	0.3	ZFX11ZB-CE50ASZ03N
		1.3	ZFX11ZB-CE50AS13N

### Gold-Tite Hexalobular Screw

### Hexalobular Try-In Screw



Part No.		Part No.	
ZFX09-ZB-CE-HLGTN	20 Ncm	ZFX09-ZB-CE-HLTISN	Hand-tighten

## ZimVie Dental Eztetic and TSX Restorative Platform Connection

### TiBase Engaging

incl. Hexalobular Retaining Screw and Gluing Screw



Post Height <sup>1</sup> (mm)	Restorative Platform (mm)	Collar Height (mm)	Part No.
7	2.9	1.0	ZFX11ZB-EZ29AS10E
		1.7	ZFX11ZB-EZ29AS17E
		2.6	ZFX11ZB-EZ29AS26E

### Hexalobular Retaining Screw



Part No.	
ZFX09-ZB-EZ-HLRS	30 Ncm

### TiBase Non-Engaging

incl. Hexalobular Retaining Screw and Gluing Screw




Post Height <sup>1</sup> (mm)	Restorative Platform (mm)	Collar Height (mm)	Part No.
7	2.9	1.0	ZFX11ZB-EZ29AS10N
		1.7	ZFX11ZB-EZ29AS17N

### Hexalobular Retaining Screw



Part No.	
ZFX09-ZB-EZ-HLRS	30 Ncm

<sup>1</sup> Various groove heights at 4.5, 5.7 mm



# NEW: GenTek TiBases

## ZimVie Dental TSX, Tapered Screw-Vent, & Trabecular Metal Restorative Platform Connections

### TiBase Engaging

*incl. Hexalobular Retaining Screw and Gluing Screw*



Post Height <sup>1</sup> (mm)	Restorative Platform (mm)	Collar Height (mm)	Part No.	
7	3.5	0.3	ZFX11ZB-TSV35AS03E	
		1.3	ZFX11ZB-TSV35AS13E	
		2.6	ZFX11ZB-TSV35AS26E	
	4.5	4.5	0.3	ZFX11ZB-TSV45AS03E
			1.3	ZFX11ZB-TSV45AS13E
			2.6	ZFX11ZB-TSV45AS26E
	5.7 <sup>1</sup>	5.7 <sup>1</sup>	0.3	ZFX11ZB-TSV57AS03E
			1.3	ZFX11ZB-TSV57AS13E
			2.6	ZFX11ZB-TSV57AS26E

### Hexalobular Retaining Screw



Part No.	Recommended Torque
ZFX09-ZB-TSV-HLRS	30 Ncm

### TiBase Non-Engaging

*incl. Hexalobular Retaining Screw and Gluing Screw*



Post Height <sup>2</sup> (mm)	Restorative Platform (mm)	Collar Height (mm)	Part No.
7	3.5	0.3	ZFX11ZB-TSV35AS03N
		1.3	ZFX11ZB-TSV35AS13N
	4.5	0.3	ZFX11ZB-TSV45AS03N
		1.3	ZFX11ZB-TSV45AS13N
	5.7 <sup>1</sup>	0.3	ZFX11ZB-TSV57AS03N
		1.3	ZFX11ZB-TSV57AS13N

### Hexalobular Retaining Screw



Part No.	Recommended Torque
ZFX09-ZB-TSV-HLRS	30 Ncm

## Tools for New GenTek TiBases

One driver for compatible ZimVie Dental Connections

### Hexalobular Screwdriver for Hexalobular Retaining Screw



Length (mm)	Part No.
21	ZFX02HLD21
28	ZFX02HLD28
36	ZFX02HLD36

### Low Torque Indicating Ratchet Wrench



Description	Part No.
Low Torque Indicating Ratchet Wrench	L-TIRW
L-TIRW Standard ISO 1797 Adapter	C9980

### Certain QuickSeat® Activator Tool



Part No.
IQSA01

### Low Torque Indicating Ratchet Wrench Kit<sup>3</sup> incl. Adapter, Tray, and Driver Tips



Description	Part No.
<ul style="list-style-type: none"> <li>L-TIRW: Low Torque Indicating Ratchet Wrench</li> <li>C9980: L-TIRW Standard ISO 1797 Adapter<sup>4</sup><sup>5</sup></li> <li>RTI2035TR: Torque Indicating Ratchet Wrench Tray</li> <li>RASH3N: Narrow Right Angle Large Hexed Driver Tip, 24mm</li> <li>RASQ3N: Narrow Right Angle Square Driver Tip, 24mm</li> </ul>	L-TIRWK

<sup>1</sup> Not compatible with TSX Restorative Platform

<sup>2</sup> Various groove heights at 4.5, 5.7mm

<sup>3</sup> Components also sold separately

<sup>4</sup> Short ISO 1797 adapter also available (C9981)

<sup>5</sup> C9980 can be ordered separately

# GenTek Pre-Abutment Blanks

## ZimVie Dental Certain Restorative Platform Connection

### Titanium Pre-milled Abutment Blank *incl. Gold-Tite Screw*



Restorative Platform (mm)	Part No.
3.4	ZFX18-ZB-CE-34
4.1, 5.0, 6.0	ZFX18-ZB-CE-4160

### Gold-Tite Screw



Part No.	Recommended Torque
ZFX09-ZB-CE-HGTS	20 Ncm

### Try-In Screw



Part No.	Recommended Torque
ZFX09-ZB-CE-HTIS	Hand-tighten

### Hex Driver, Driver Tip *for Gold-Tite and Try-In Screws*



Description	Part No.
Narrow Posterior Hex Driver, 17 mmL	PHD02N
Narrow Standard Hex Driver, 24 mmL	PHD03N
Narrow Right Angle Hex Driver Tip, 24 mmL	RASH3N
Narrow Right Angle Hex Driver Tip, 30 mmL	RASH8N

### Low Torque Indicating Ratchet Wrench Kit<sup>1</sup> *incl. Adapter, Tray, and Driver Tips*



Description	Part No.
<ul style="list-style-type: none"> <li>L-TIRW: Low Torque Indicating Ratchet Wrench</li> <li>C9980: L-TIRW Standard ISO 1797 Adapter<sup>2/3</sup></li> <li>RTI2035TR: Torque Indicating Ratchet Wrench Tray</li> <li>RASH3N: Narrow Right Angle Large Hexed Driver Tip, 24mm</li> <li>RASQ3N: Narrow Right Angle Square Driver Tip, 24mm</li> </ul>	L-TIRWK

### Certain QuickSeat Activator Tool



Part No.
IQSA01



## ZimVie Dental External Hex Restorative Platform Connection

### Titanium Pre-milled Abutment Blank *incl. Gold-Tite Screw*



Restorative Platform (mm)	Part No.
3.4	ZFX18-ZB-EX-34
4.1, 5.0, 6.0	ZFX18-ZB-EX-4160

### Gold-Tite Screw



Part No.	Recommended Torque
ZFX09-ZB-EX-GTS	32 Ncm

### Try-In Screw



Part No.	Recommended Torque
ZFX09-ZB-EX-TIS	Hand-tighten

### Square Driver, Driver Tip *for Gold-Tite Screw*



Description	Part No.
Narrow Posterior Square Driver, 17 mmL	PSQD0N
Narrow Standard Slotted Square Driver, 24 mmL	PSQD1N
Narrow Right Angle Square Driver Tip, 24 mmL	RASQ3N
Narrow Right Angle Square Driver Tip, 30 mmL	RASH8N

### Hex Driver, Driver Tip *for Try-In Screw*



Only for External Hex Try-In Screw

Description	Part No.
Narrow Posterior Hex Driver, 17 mmL	PHD02N
Narrow Standard Hex Driver, 24 mmL	PHD03N
Narrow Right Angle Hex Driver Tip, 24 mmL	RASH3N
Narrow Right Angle Hex Driver Tip, 30 mmL	RASH8N

### Low Torque Indicating Ratchet Wrench Kit<sup>1</sup> *incl. Adapter, Tray, and Driver Tips*




Description	Part No.
<ul style="list-style-type: none"> <li>L-TIRW: Low Torque Indicating Ratchet Wrench</li> <li>C9980: L-TIRW Standard ISO 1797 Adapter<sup>2,3</sup></li> <li>RTI2035TR: Torque Indicating Ratchet Wrench Tray</li> <li>RASH3N: Narrow Right Angle Large Hexed Driver Tip, 24 mm</li> <li>RASQ3N: Narrow Right Angle Square Driver Tip, 24 mm</li> </ul>	L-TIRWK

<sup>1</sup>Components also sold separately

<sup>2</sup>Short ISO 1797 adapter also available (C9981)

<sup>3</sup>C9980 can be ordered separately



# GenTek Pre-Abutment Blanks

## ZimVie Dental Eztetic and TSX Restorative Platform Connection

### Titanium Pre-milled Abutment Blank *incl. Retaining Screw*



Restorative Platform (mm)	Part No.	Part No.	
2.9	ZFX18-ZB-EZ-29	ZFX09-ZB-EZ-RS	30 Ncm

### Retaining Screw




### Hex Driver, Latch-Lock Hex Driver *for Retaining Screws & Intraoral Scanbodies*



Description	Part No.
Hex Driver, Short, with GemLock® Retention 1.25 mmD, 22 mmL	HXGR1.25 <sup>1 2 3</sup>
Hex Driver, Long, with GemLock Retention 1.25 mmD, 30 mmL	HXLGR1.25 <sup>1 2 3</sup>
Latch-Lock Hex Driver for Surgical and Prosthetic Screws, Short 1.25 mmD, 23 mmL <i>(may require reduction hand piece)</i>	HX1.25D <sup>1</sup>
Latch-Lock Hex Driver for Surgical and Prosthetic Screws, Long 1.25 mmD, 26 mmL <i>(may require reduction hand piece)</i>	HXL1.25D <sup>1</sup>

### Removal Tool



Description	Part No.
Removal Tool for Abutments	CLRT2
Removal Tool for Narrow Platform for Abutments	ZFX02-RTNP

### Torque Wrench Tool

*for Abutment Screws & Intraoral Scanbodies*



Description	Part No.
Torque Wrench, Restorative (adjustable torque range 10 – 35 Ncm)	TWR
Torque Wrench Driver, Short	TW1.25 <sup>1</sup>
Torque Wrench Driver, Long	TW1.25L <sup>1 2</sup>

## ZimVie Dental TSX, Tapered Screw-Vent, & Trabecular Metal Restorative Platform Connection

### Titanium Pre-milled Abutment Blank *incl. Retaining Screw*



### Retaining Screw




Restorative Platform (mm)	Part No.	Part No.	
3.5	ZFX18-ZB-TSV-35	ZFX09-ZB-TSV-RS	30 Ncm
4.5	ZFX18-ZB-TSV-45		
5.7 <sup>5</sup>	ZFX18-ZB-TSV-57		

### Hex Driver, Latch-Lock Hex Driver *for Retaining Screws & Intraoral Scanbodies*



Description	Part No.
Hex Driver, Short, with GemLock Retention 1.25 mmD, 22 mmL	HXGR1.25 <sup>1 2 4</sup>
Hex Driver, Long, with GemLock Retention 1.25 mmD, 30 mmL	HXLGR1.25 <sup>1 2 4</sup>
Latch-Lock Hex Driver for Surgical and Prosthetic Screws, Short 1.25 mmD, 23 mmL <i>(may require reduction hand piece)</i>	HX1.25D <sup>1</sup>
Latch-Lock Hex Driver for Surgical and Prosthetic Screws, Long 1.25 mmD, 26 mmL <i>(may require reduction hand piece)</i>	HXL1.25D <sup>1</sup>

### Removal Tool



Description	Part No.
Removal Tool for Abutments	TLRT2
Removal Tool Standard Platform for Abutments	ZFX02-RTSP

### Torque Wrench Tool

*for Abutment Screws & Intraoral Scanbodies*



Description	Part No.
Torque Wrench, Restorative (adjustable torque range 10 – 35 Ncm)	TWR
Torque Wrench Driver, Short	TW1.25 <sup>1</sup>
Torque Wrench Driver, Long	TW1.25L <sup>1 2</sup>

<sup>1</sup> Compatible with TiBases

<sup>2</sup> Compatible with IOS Scanbodies

<sup>3</sup> Compatible with Removal Tool Narrow Platform

<sup>4</sup> Compatible with Removal Tool Standard Platform

<sup>5</sup> Not compatible with TSX Restorative Platform

# GenTek TiBases Milling Block Compatible

## ZimVie Dental Certain Restorative Platform Connection

### Single-Unit Restorative Components

**TiBase Engaging<sup>1</sup>**  
incl. Gold-Tite and Gluing Screws



Height (mm)	Restorative Platform (mm)	Part No.
4.7	3.4	ZFX11-ZB-CE-3447-ES
	4.1	ZFX11-ZB-CE-4147-EL
	5.0	ZFX11-ZB-CE-5047-EL

**Gold-Tite Screw**

**Try-In Screw**



Part No.		Part No.	
ZFX09-ZB-CE-HGTS	20 Ncm	ZFX09-ZB-CE-HTIS	Hand-tighten

### Multi-Unit Restorative Components

**TiBase Non-Engaging<sup>1</sup>**  
incl. Gold-Tite and Gluing Screws



Height (mm)	Restorative Platform (mm)	Part No.
4.7	3.4	ZFX11-ZB-CE-3447-NES
	4.1	ZFX11-ZB-CE-4147-NEL
	5.0	ZFX11-ZB-CE-5047-NEL

**Gold-Tite Screw**

**Try-In Screw**



Part No.		Part No.	
ZFX09-ZB-CE-LGTS	20 Ncm	ZFX09-ZB-CE-NHTIS	Hand-tighten

**Hex Driver, Driver Tip**  
for Gold-Tite and Try-In Screws



Description	Part No.
Narrow Posterior Hex Driver, 17 mmL	PHD02N
Narrow Standard Hex Driver, 24 mmL	PHD03N
Narrow Right Angle Hex Driver Tip, 24 mmL	RASH3N
Narrow Right Angle Hex Driver Tip, 30 mmL	RASH8N

**Low Torque Indicating Ratchet Wrench Kit<sup>2</sup>**  
incl. Adapter, Tray, and Driver Tips



**Certain QuickSeat<sup>®</sup> Activator Tool**



Description	Part No.	Part No.
<ul style="list-style-type: none"> <li>L-TIRW: Low Torque Indicating Ratchet Wrench</li> <li>C9980: L-TIRW Standard ISO 1797 Adapter<sup>3,4</sup></li> <li>RTI2035TR: Torque Indicating Ratchet Wrench Tray</li> <li>RASH3N: Narrow Right Angle Large Hexed Driver Tip, 24 mm</li> <li>RASQ3N: Narrow Right Angle Square Driver Tip, 24 mm</li> </ul>	L-TIRWK	IQSA01


<sup>1</sup>4 mm groove indicating minimal allowed post height

<sup>2</sup>Not applicable for the GenTek Angulated Screw Channel TiBase

<sup>3</sup>Components also sold separately

<sup>4</sup>Short ISO 1797 adapter also available (C9981)

<sup>5</sup>C9980 can be ordered separately



# GenTek TiBases Milling Block Compatible

## ZimVie Dental External Hex Restorative Platform Connection

### Single-Unit Restorative Components

**TiBase Engaging<sup>1</sup>**  
*incl. Gold-Tite and Gluing Screws*



Height (mm)	Restorative Platform (mm)	Part No.
4.7	3.4	ZFX11-ZB-EX-3447-EL
	4.1	ZFX11-ZB-EX-4147-EL
	5.0	ZFX11-ZB-EX-5047-EL

**Gold-Tite Screw**



Recommended Torque

Part No.	
ZFX09-ZB-EX-GTS	32 Ncm

**Try-In Screw**



Recommended Torque

Part No.	
ZFX09-ZB-EX-TIS	Hand-tighten

### Multi-Unit Restorative Components

**TiBase Non-Engaging<sup>1</sup>**  
*incl. Gold-Tite and Gluing Screws*



Height (mm)	Restorative Platform (mm)	Part No.
4.7	3.4	ZFX11-ZB-EX-3447-NEL
	4.1	ZFX11-ZB-EX-4147-NEL
	5.0	ZFX11-ZB-EX-5047-NEL

**Gold-Tite Screw**



Recommended Torque

Part No.	
ZFX09-ZB-EX-GTS	32 Ncm

**Try-In Screw**



Recommended Torque

Part No.	
ZFX09-ZB-EX-TIS	Hand-tighten

**Square Driver, Driver Tip**  
for Gold-Tite Screw



**Hex Driver, Driver Tip**  
for Try-In Screw



Only for  
External Hex  
Try-In Screw

Description	Part No.
Narrow Posterior Square Driver, 17 mmL	PSQD0N
Narrow Standard Slotted Square Driver, 24 mmL	PSQD1N
Narrow Right Angle Square Driver Tip, 24 mmL	RASQ3N
Narrow Right Angle Square Driver Tip, 30 mmL	RASH8N

Description	Part No.
Narrow Posterior Hex Driver, 17 mmL	PHD02N
Narrow Standard Hex Driver, 24 mmL	PHD03N
Narrow Right Angle Hex Driver Tip, 24 mmL	RASH3N
Narrow Right Angle Hex Driver Tip, 30 mmL	RASH8N

**Low Torque Indicating Ratchet Wrench Kit<sup>2</sup>**  
incl. Adapter, Tray, and Driver Tips



Description	Part No.
<ul style="list-style-type: none"> <li>L-TIRW: Low Torque Indicating Ratchet Wrench</li> <li>C9980: L-TIRW Standard ISO 1797 Adapter<sup>3,4</sup></li> <li>RTI2035TR: Torque Indicating Ratchet Wrench Tray</li> <li>RASH3N: Narrow Right Angle Large Hexed Driver Tip, 24 mm</li> <li>RASQ3N: Narrow Right Angle Square Driver Tip, 24 mm</li> </ul>	L-TIRWK

<sup>1</sup>4 mm groove indicating minimal allowed post height

<sup>2</sup>Not applicable for the GenTek Angulated Screw Channel TiBase

<sup>3</sup>Components also sold separately

<sup>4</sup>Short ISO 1797 adapter also available (C9981)

<sup>5</sup>C9980 can be ordered separately

# GenTek TiBases Milling Block Compatible

## ZimVie Dental Eztetic and TSX Restorative Platform Connection

### Single-Unit Restorative Components

**TiBase Engaging<sup>1</sup>**  
incl. Retaining and Gluing Screws



Height (mm)	Restorative Platform (mm)	Part No.
4.7	2.9	ZFX11-ZB-EZ-2947-ES

### Retaining Screw



Part No.	Recommended Torque
ZFX09-ZB-EZ-RS	30 Ncm

### Multi-Unit Restorative Components

**TiBase Non-Engaging<sup>1</sup>**  
incl. Retaining and Gluing Screws



Height (mm)	Restorative Platform (mm)	Part No.
4.7	2.9	ZFX11-ZB-EZ-2947-NES

### Retaining Screw



Part No.	Recommended Torque
ZFX09-ZB-EZ-RS	30 Ncm



**Hex Driver, Latch-Lock Hex Driver**  
for Retaining Screws



Description	Part No.
Hex Driver, Short, with GemLock Retention 1.25 mmD, 22 mmL	HXGR1.25 <sup>2 3 4</sup>
Hex Driver, Long, with GemLock Retention 1.25 mmD, 30 mmL	HXLGR1.25 <sup>2 3 4</sup>
Latch-Lock Hex Driver for Surgical and Prosthetic Screws, Short 1.25 mmD, 23 mmL (may require reduction hand piece)	HX1.25D <sup>2</sup>
Latch-Lock Hex Driver for Surgical and Prosthetic Screws, Long 1.25 mmD, 26 mmL (may require reduction hand piece)	HXL1.25D <sup>2</sup>

**Removal Tool**



Description	Part No.
Removal Tool for Abutments	CLRT2
Removal Tool for Narrow Platform for Abutments	ZFX02-RTNP


**Torque Wrench Tool**  
for Abutment Screws



Description	Part No.
Torque Wrench, Restorative (adjustable torque range 10 – 35 Ncm)	TWR
Torque Wrench Driver, Short	TW1.25 <sup>2</sup>
Torque Wrench Driver, Long	TW1.25L <sup>2 3</sup>

<sup>1</sup>4 mm groove indicating minimal allowed post height  
Not applicable for the GenTek Angulated Screw Channel TiBase  
<sup>2</sup>Compatible with TiBases

<sup>3</sup>Compatible with IOS Scanbodies  
<sup>4</sup>Compatible with Removal Tool Narrow Platform



# GenTek TiBases Milling Block Compatible

## ZimVie Dental TSX, Tapered Screw-Vent, & Trabecular Metal Restorative Platform Connections

### Single-Unit Restorative Components

#### TiBase Engaging<sup>1</sup> incl. Retaining and Gluing Screws



Height (mm)	Restorative Platform (mm)	Part No.
4.7	3.5	ZFX11-ZB-TSV-3547-EL
	4.5	ZFX11-ZB-TSV-4547-EL
	5.7 <sup>3</sup>	ZFX11-ZB-TSV-5747-EL

#### Retaining Screw



Part No.	
ZFX09-ZB-TSV-RS	30 Ncm

### Multi-Unit Restorative Components

#### TiBase Non-Engaging<sup>1</sup> incl. Retaining and Gluing Screws



Height (mm)	Restorative Platform (mm)	Part No.
4.7	3.5	ZFX11-ZB-TSV-3547-NEL
	4.5	ZFX11-ZB-TSV-4547-NEL
	5.7 <sup>3</sup>	ZFX11-ZB-TSV-5747-NEL

#### Retaining Screw




Part No.	
ZFX09-ZB-TSV-RS	30 Ncm

**Hex Driver, Latch-Lock Hex Driver**  
for Retaining Screws



Description	Part No.
Hex Driver, Short, with GemLock Retention 1.25 mmD, 22 mmL	HXGR1.25 <sup>2 4 5</sup>
Hex Driver, Long, with GemLock Retention 1.25 mmD, 30 mmL	HXLGR1.25 <sup>2 4 5</sup>
Latch-Lock Hex Driver for Surgical and Prosthetic Screws, Short 1.25 mmD, 23 mmL (may require reduction hand piece)	HX1.25D <sup>2</sup>
Latch-Lock Hex Driver for Surgical and Prosthetic Screws, Long 1.25 mmD, 26 mmL (may require reduction hand piece)	HXL1.25D <sup>2</sup>

**Removal Tool**



Description	Part No.
Removal Tool for Abutments	TLRT2
Removal Tool Standard Platform for Abutments	ZFX02-RTSP

**Torque Wrench Tool**  
for Abutment Screws



Description	Part No.
Torque Wrench, Restorative (adjustable torque range 10 – 35 Ncm)	TWR
Torque Wrench Driver, Short	TW1.25 <sup>2</sup>
Torque Wrench Driver, Long	TW1.25L <sup>2 4</sup>

<sup>1</sup>4 mm groove indicating minimal allowed post height

<sup>5</sup>Not applicable for the GenTek Angulated Screw Channel TiBase

<sup>2</sup>Compatible with TiBases

<sup>3</sup>Not compatible with TSX Restorative Platform

<sup>4</sup>Compatible with IOS Scanbodies

<sup>5</sup>Compatible with Removal Tool Standard Platform

# GenTek TiBases ZimVie Multi-Unit Abutment Compatible

## ZimVie Dental Low Profile Abutment

### Multi-Unit Restorative Components

#### TiBase Non-Engaging<sup>1</sup> incl. Gold-Tite and Gluing Screws



Height (mm)	Restorative Platform (mm)	Part No.
5.0	4.8	ZFX11-ZB-LP-48

#### Gold-Tite Screw

#### Try-in Screw



Part No.		Recommended Torque	Part No.		Recommended Torque
ZFX09-ZB-LP-GTS		10 Ncm	ZFX09-ZB-LP-RS		Hand-tighten

#### Hex Driver, Driver Tip for Gold-Tite Screws and Try-in Screws



Description	Part No.
Narrow Posterior Hex Driver, 17 mmL	PHD02N
Narrow Standard Hex Driver, 24 mmL	PHD03N
Narrow Right Angle Hex Driver Tip, 24 mmL	RASH3N
Narrow Right Angle Hex Driver Tip, 30 mmL	RASH8N

#### Low Torque Indicating Ratchet Wrench Kit<sup>2</sup> incl. Adapter, Tray, and Driver Tips



Description	Part No.
<ul style="list-style-type: none"> <li>L-TIRW: Low Torque Indicating Ratchet Wrench</li> <li>C9980: L-TIRW Standard ISO 1797 Adapter<sup>3,4</sup></li> <li>RTI2035TR: Torque Indicating Ratchet Wrench Tray</li> <li>RASH3N: Narrow Right Angle Large Hexed Driver Tip, 24 mm</li> <li>RASQ3N: Narrow Right Angle Square Driver Tip, 24 mm</li> </ul>	L-TIRWK

## ZimVie Dental Tapered Abutment

### Multi-Unit Restorative Components

#### TiBase Non-Engaging<sup>5</sup> incl. Retaining and Gluing Screws



Height (mm)	Restorative Platform (mm)	Part No.
3.9	4.5	ZFX11-ZB-TA-45

#### Retaining Screw



Part No.	Recommended Torque
ZFX09-ZB-TA-RS	20 Ncm

#### Hex Driver, Latch-Lock Hex Driver for Retaining Screws



Description	Part No.
Hex Driver, Short, with GemLock Retention 1.25 mmD, 22 mmL	HXGR1.25
Hex Driver, Long, with GemLock Retention 1.25 mmD, 30 mmL	HXLGR1.25
Latch-Lock Hex Driver for Surgical and Prosthetic Screws, Short 1.25 mmD, 23 mmL (may require reduction hand piece)	HX1.25D
Latch-Lock Hex Driver for Surgical and Prosthetic Screws, Long 1.25 mmD, 26 mmL (may require reduction hand piece)	HXL1.25D

#### Torque Wrench Tool for Abutment Screws



Description	Part No.
Torque Wrench, Restorative (adjustable torque range 10 – 35 Ncm)	TWR
Torque Wrench Driver, Short	TW1.25 <sup>6</sup>
Torque Wrench Driver, Long	TW1.25L <sup>6,7</sup>

<sup>1</sup> Compatible with non-hexed low profile abutment only  
<sup>2</sup> Not applicable for the GenTek Angulated Screw Channel TiBase

<sup>3</sup> Components also sold separately  
<sup>4</sup> Short ISO 1797 adapter also available (C9981)  
<sup>5</sup> C9980 can be ordered separately

<sup>6</sup> Not applicable for the GenTek Angulated Screw Channel TiBase  
<sup>7</sup> Compatible with TiBases  
<sup>8</sup> Compatible with IOS Scanbodies



## ZimVie Dental Certain Restorative Platform Connection

### Desktop Scanbody incl. Screw



Restorative Platform (mm)	Part No.	
3.4	ZFX05D-ZB-CE-34	Hand-tighten <sup>1</sup>
4.1, 5.0, 6.0	ZFX05D-ZB-CE-4160	
5.0	ZFX05D-ZB-CE-50 <sup>2</sup>	
6.0	ZFX05D-ZB-CE-60 <sup>2</sup>	

### Torque Wrench for Desktop Scanbody



Part No.
ZFX02002067

### Intraoral Scanbodies incl. Screws



Height (mm)	Restorative Platform (mm)	Part No.	
4.0 & 7.0	3.4	ZFX05I-ZB-CE-34	Hand-tighten <sup>1</sup>
	4.1, 5.0, 6.0	ZFX05I-ZB-CE-4160	

### Hex Driver for Intraoral Scanbodies



Description	Part No.
Hex Driver, Long, with GemLock Retention 1.25 mmD, 30 mmL	HXLGR1.25

### Certain QuickSeat® Activator Tool



Part No.
IQSA01

## ZimVie Dental External Hex Restorative Platform Connection

### Desktop Scanbody incl. Screw



↑  
Recommended  
Torque

Restorative Platform (mm)	Part No.	Hand-tighten <sup>1</sup>
3.4	ZFX05D-ZB-EX-34	
4.1, 5.0, 6.0	ZFX05D-ZB-EX-4160	
5.0	ZFX05D-ZB-EX-50 <sup>2</sup>	
6.0	ZFX05D-ZB-EX-60 <sup>2</sup>	

### Torque Wrench for Desktop Scanbody



Part No.
ZFX02002067

### Intraoral Scanbodies incl. Screws



↑  
Recommended  
Torque


Height (mm)	Restorative Platform (mm)	Part No.	Hand-tighten <sup>1</sup>
4.0 & 7.0	3.4	ZFX05I-ZB-EX-34	
	4.1, 5.0, 6.0	ZFX05I-ZB-EX-4160	

### Hex Driver for Intraoral Scanbodies



Description	Part No.
Hex Driver, Long, with GemLock Retention 1.25 mmD, 30 mmL	HXLGR1.25

<sup>1</sup>Minimum of 5 Ncm; Maximum of 15 Ncm.  
Refer to Instruction Manual for more information.  
<sup>2</sup>Only for use with Zfx Evolution scanners



## ZimVie Dental Eztetic and TSX Restorative Platform Connection

### Desktop Scanbody incl. Screw



↑  
Recommended  
Torque


Restorative Platform (mm)	Part No.	
2.9	ZFX05D-ZB-EZ-29	Hand-tighten <sup>1</sup>

### Torque Wrench for Desktop Scanbody



Part No.
ZFX02002067

### Intraoral Scanbodies incl. Screws



↑  
Recommended  
Torque

Height (mm)	Restorative Platform (mm)	Part No.	
4.0 & 7.0	2.9	ZFX05I-ZB-EZ-29	Hand-tighten <sup>1</sup>

### Hex Driver, Latch-Lock Hex Driver for Intraoral Scanbodies



Description	Part No.
Hex Driver, Short, with GemLock Retention 1.25 mmD, 22 mmL	HXGR1.25 <sup>2 3 4</sup>
Hex Driver, Long, with GemLock Retention 1.25 mmD, 30 mmL	HXLGR1.25 <sup>2 3 4</sup>
Latch-Lock Hex Driver for Surgical and Prosthetic Screws, Short 1.25 mmD, 23 mmL (may require reduction hand piece)	HX1.25D <sup>2</sup>
Latch-Lock Hex Driver for Surgical and Prosthetic Screws, Long 1.25 mmD, 26 mmL (may require reduction hand piece)	HXL1.25D <sup>2</sup>



## ZimVie Dental TSX, Tapered Screw-Vent, & Trabecular Metal Restorative Platform Connection

### Desktop Scanbody incl. Screw



↑  
Recommended  
Torque

Restorative Platform (mm)	Part No.	Hand-tighten <sup>1</sup>
3.5	ZFX05D-ZB-TSV-35	
4.5	ZFX05D-ZB-TSV-45	
5.7 <sup>6</sup>	ZFX05D-ZB-TSV-57	

### Torque Wrench for Desktop Scanbody



Part No.
ZFX02002067

### Intraoral Scanbodies incl. Screws



↑  
Recommended  
Torque

Height (mm)	Restorative Platform (mm)	Part No.	Hand-tighten <sup>1</sup>
4.0 & 7.0	3.5	ZFX05I-ZB-TSV-35	
	4.5	ZFX05I-ZB-TSV-45	
	5.7 <sup>6</sup>	ZFX05I-ZB-TSV-57	

### Hex Driver, Latch-Lock Hex Driver for Intraoral Scanbodies



Description	Part No.
Hex Driver, Short, with GemLock Retention 1.25 mmD, 22 mmL	HXGR1.25 <sup>2 3 5</sup>
Hex Driver, Long, with GemLock Retention 1.25 mmD, 30 mmL	HXLGR1.25 <sup>2 3 5</sup>
Latch-Lock Hex Driver for Surgical and Prosthetic Screws, Short 1.25 mmD, 23 mmL (may require reduction hand piece)	HX1.25D <sup>2</sup>
Latch-Lock Hex Driver for Surgical and Prosthetic Screws, Long 1.25 mmD, 26 mmL (may require reduction hand piece)	HXL1.25D <sup>2</sup>

<sup>1</sup> Minimum of 5 Ncm; Maximum of 15 Ncm. Refer to Instruction Manual for more information.

<sup>2</sup> Compatible with TiBases

<sup>3</sup> Compatible with IOS Scanbodies

<sup>4</sup> Compatible with Removal Tool Narrow Platform


<sup>5</sup> Compatible with Removal Tool Standard Platform

<sup>6</sup> Not compatible with TSX Restorative Platform



**ZimVie Dental Low Profile Abutment**

**Desktop Scanbody**  
*incl. Screw*



↑  
Recommended  
Torque

Restorative Platform (mm)	Part No.		
4.8	ZFX05D-ZB-LP-48		Hand-tighten <sup>1</sup>

**Torque Wrench**  
*for Desktop Scanbody*



Part No.
ZFX02002067

**Intraoral Scanbodies Engaging**  
*incl. Screws*



↑  
Recommended  
Torque

Height (mm)	Restorative Platform (mm)	Part No.	
4.0	4.8	ZFX05I-ZB-LP-48	Hand-tighten <sup>1</sup>

**Hex Driver, Latch-Lock Hex Driver**  
*for Intraoral Scanbodies*



Description	Part No.
Hex Driver, Short, with GemLock Retention 1.25 mmD, 22 mmL	HXGR1.25
Hex Driver, Long, with GemLock Retention 1.25 mmD, 30 mmL	HXLGR1.25
Latch-Lock Hex Driver for Surgical and Prosthetic Screws, Short, 1.25 mmD, 23 mmL <i>(may require reduction hand piece)</i>	HX1.25D
Latch-Lock Hex Driver for Surgical and Prosthetic Screws, Long, 1.25 mmD, 26 mmL <i>(may require reduction hand piece)</i>	HXL1.25D

**Intraoral Scanbodies Non-Engaging**  
*incl. Screws*




↑  
Recommended  
Torque

Height (mm)	Restorative Platform (mm)	Part No.	
4.0	4.8	ZFX05I-ZB-LP-48-NE	Hand-tighten <sup>1</sup>

## ZimVie Dental Tapered Abutment

### Desktop Scanbody incl. Screw



↑  
Recommended  
Torque


Restorative Platform (mm)	Part No.	Hand-tighten <sup>1</sup>
4.5	ZFX05D-ZB-TA-45	

### Torque Wrench for Desktop Scanbody



Part No.
ZFX02002067

### Intraoral Scanbodies incl. Screws



DOUBLE PACK

↑  
Recommended  
Torque


Height (mm)	Restorative Platform (mm)	Part No.	Hand-tighten <sup>1</sup>
4.0	4.5	ZFX05I-ZB-TA-45	

### Hex Driver, Latch-Lock Hex Driver for Intraoral Scanbodies



Description	Part No.
Hex Driver, Short, with GemLock Retention 1.25 mmD, 22 mmL	HXGR1.25
Hex Driver, Long, with GemLock Retention 1.25 mmD, 30 mmL	HXLGR1.25
Latch-Lock Hex Driver for Surgical and Prosthetic Screws, Short, 1.25 mmD, 23 mmL (may require reduction hand piece)	HX1.25D
Latch-Lock Hex Driver for Surgical and Prosthetic Screws, Long, 1.25 mmD, 26 mmL (may require reduction hand piece)	HXL1.25D

<sup>1</sup> Minimum of 5 Ncm; Maximum of 15 Ncm.  
Refer to Instruction Manual for more information.



# GenTek Digital Scan Analogs

## ZimVie Dental Certain Restorative Platform Connection

### Digital Scan Analog<sup>1</sup>



Restorative Platform (mm)	Part No.
3.4	ZFX15-ZB-CE-34
4.1	ZFX15-ZB-CE-41
5.0	ZFX15-ZB-CE-50
6.0	ZFX15-ZB-CE-60

### Placement Tool for Digital Scan Analogs



Part No.
ZFX02-M16

## ZimVie Dental External Hex Restorative Platform Connection

### Digital Scan Analog<sup>1</sup>



Restorative Platform (mm)	Part No.
3.4	ZFX15-ZB-EX-34
4.1	ZFX15-ZB-EX-41
5.0	ZFX15-ZB-EX-50
6.0	ZFX15-ZB-EX-60

### Placement Tool for Digital Scan Analogs



Part No.
ZFX02-M20

## ZimVie Dental Low Profile Abutment

### Digital Scan Analog<sup>1</sup>



Restorative Platform (mm)	Part No.
4.8	ZFX15-ZB-LP-48

### Placement Tool for Digital Scan Analogs



Part No.
ZFX02-M14

## ZimVie Dental TSX, Tapered Screw-Vent, & Trabecular Metal Restorative Platform Connection

### Digital Scan Analog<sup>1</sup>



### Placement Tool for Digital Scan Analogs



Restorative Platform (mm)	Part No.
3.5	ZFX15-ZB-TSV-35
4.5	ZFX15-ZB-TSV-45
5.7	ZFX15-ZB-TSV-57 <sup>2</sup>

Part No.
ZFX02-M18

## ZimVie Dental Eztetic and TSX Restorative Platform Connection

### Digital Scan Analog<sup>1</sup>



### Placement Tool for Digital Scan Analogs



Restorative Platform (mm)	Part No.
2.9	ZFX15-ZB-EZ-29

Part No.
ZFX02-M16

## ZimVie Dental Tapered Abutment

### Digital Scan Analog<sup>1</sup>



### Placement Tool for Digital Scan Analogs



Restorative Platform (mm)	Part No.
4.5	ZFX15-ZB-TA-45

Part No.
ZFX02-M18

<sup>1</sup>Compatible with existing Certain / External Hex / Eztetic / TSV/TM / Low Profile / Tapered Abutment Impression Copings

<sup>2</sup>Not compatible with TSX Restorative Platform





To order genuine ZimVie Dental Components  
call your local representative or send an email to  
**ZV.UKorders@ZimVie.com**

**ZimVie Dental Global Headquarters**

4555 Riverside Drive  
Palm Beach Gardens, FL 33410  
Phone: +1-561-776-6700  
Fax: +1-561-776-1272  
dentalCS@ZimVie.com  
www.ZimVie.com/dental

**Biomet 3i (UK & Ireland) Ltd.**

Reading Business Centre  
Suite 807, 8th Floor Fountain House  
2 Queens Walk, Reading, Berks  
RG1 7QF, United Kingdom  
Telephone (UK): 0800 652 1233  
Telephone (Ireland): 1 800 552 752  
ZV.UKorders@ZimVie.com



Unless otherwise indicated, as referenced herein, all trademarks and intellectual property rights are the property of ZimVie Inc. or an affiliate; and all products are manufactured by one or more of the dental subsidiaries of ZimVie Inc. (Biomet 3i, LLC, Zimmer Dental, Inc., etc.) and marketed and distributed by ZimVie Dental and its authorized marketing partners. For additional product information, please refer to the individual product labeling or instructions for use. Product clearance and availability may be limited to certain countries/regions. This material is intended for clinicians only and does not comprise medical advice or recommendations. Distribution to any other recipient is prohibited. This material may not be copied or reprinted without the express written consent of ZimVie. TiBases (Class IIb), Scanbodies (Class I), Blanks (Class IIb), Tools (Class I), Screws (Class IIb) and Analogs (non medical devices). ZV0121 REV B 05/23 ©2023 ZimVie. All rights reserved.

